



	<b>VOCO</b> Anton-Flettner-Straße 1-3, 27472 Cuxhaven, Germany Tel: +49 (0)4721-719-0 Web: www.voco.com Email: c.luo@voco.de Dr. Luo Chuan Hao	<b>Heraeus Dental</b> 10 Ang Mo Kio #06-05 Techpoint Singapore 569059 Tel: +65 65717569 Web: www.heraeus-dental.com Email: dave.tan@heraeus.com Dave Tan	<b>PROMEDICA Dental Material GmbH</b> Domagkstr. 31-24537 Neumünster Tel: (+49) 43 21 – 5 41 73 Web: www.promedica.de Email: info@promedica.de Mr. Peer Duwe (Export Director)	
<b>Contact Person :</b>				
<b>BONDING</b>	<b>Product Name :</b>	Futura Bond DC	iBond® Self Etch Bonding	Compobond DC
	<b>Generation :</b>	6th Generation	7th Generation	6th Generation
	<b>Delivery :</b>	2 Bottles Bond	4 packaging	2 Bottle System
	<b>Volume:</b>	4ml+4ml	4 ml for bottle, Singles Dose 0.15ml	4 ml Liquid 1, 4 ml Liquid 2
	<b>Indications:</b>	For all classes of cavity, especially cervical filling can eliminate post operative sensitivities.	Bonding of direct and indirect composite restorations Sealing hypersensitive areas of teeth	direct self- or dual-curing composite restorations and core build-ups
	<b>Solvent:</b>	Ethanol	Acetone	BIS-GMA, HEMA, TMPTMA, campherchinon, amines (DABE), BHT, catalysts, fluorides and ethanol
	<b>Bond Strength:</b>	High bonding strength	Dentine 38 MPa (37-45) Enamel: 32 MPa (28-35)	High
	<b>Application Time:</b>	Less than 40sec	20 sec	approx. 37 seconds
	<b>Curing Time:</b>	20sec	20 sec with LED light	with blue light for 10 seconds
<b>DIRECT RESTORATIVE MATERIALS</b>	<b>Product Name :</b>	Grandio	Charisma® Opal	Light-curing nano-ceram composite
	<b>Delivery:</b>	Syringe & capsules available	Syringe	Syringe
	<b>Volume/Weight:</b>	4g syringe and 20caps/box	4 gm	4 g / syringe
	<b>Indications:</b>	Filling for all classes of cavity	Fillings of the classes I-V according to Black Shape	Fillings of class I to V
	<b>Filler Percentage:</b>	87% High filler degree	64% filler by volume	High
	<b>Curing Time:</b>	20sec	20 sec and 40 sec depending on shades	40 s
	<b>Radiopacity:</b>	Good Shrinkage; Less than 1.6%	Yes	High